**eye
banking**

DMEK vacuum
donor punch

manuale per procedura **pre-stripping**
manual for **pre-stripping** procedure

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Description

1. Punch blade $\varnothing$ 10.0 mm with calibrated depth 100 μ
2. Support
3. Syringe tubing
4. 5 cc spring syringe

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Unique features

- Cornea centering system
- 100 μ fixed-cutting depth allows safe and quick detection of the Descemet Membrane
- Predictable and reproducible system, does not depend on the operator
- No induced cellular mortality

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Pick up the cornea from the bottle and remove the Corneal Claw (suture thread); center the lenticle on the punch support using as reference the outer holes of the suction area (Fig.0). Use the syringe to create the vacuum and secure the cornea to the support.

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
The blade (Fig.1) is calibrated for a controlled penetration of 100 microns (Fig.2), thus avoiding the perforation of the stroma that could jeopardize a proper preparation. Press the blade and push it with a force perpendicular to the support without rotational movement (Fig.3), avoiding induced torsions on the membrane.

The applied force may be equal to the force applied for a perforating keratoplasty.

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Dye the endothelium with Trypan Blue or equivalent for about 20 seconds to define the cutting area.

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
Use a straight sharp forceps J3277 to remove the portion of the membrane in excess of the punching area. Re-dye using Trypan Blue or equivalent for 20 seconds and rinse with PBS, in order to highlight the portion between the membrane and the residual stroma. Check the presence of residuals which may cause endothelial adhesions during the stripping procedure. Fill the convexity of the cornea with PBS or fluid culture for storage of corneas, to the point of creating a film of fluid on top of the corneal-scleral ring, in order to avoid any tissue damages. With the Fogla cleavage hook (item number J2891E) detect the cleavage plane at the edge of the punching area (fig.4) and start lifting the membrane and carry on the procedure along the entire circumference of the cut (fig.5). Make sure there are no radial ruptures of the membrane. In the event this occurs, leave this portion at the end of stripping, so as to leave it with the membrane still adherent to the stroma. (Fig.8.1)

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Tying and stitch forceps, with delicate 6 mm tying platform, straight. For nylon and monofilament sutures. Recommended for 8-0 to 10-0 sutures

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Delicate blunt hook for Descemet membrane cleavage. The Fogla hook, with atraumatic tips, triangular shaped, performs the cleavage at the edges of the lenticle previously secured to the suction DMEK support. Once the correct thickness of the pure Descemet membrane is achieved, the Janach forceps J2892E performs the detachment by means of the peeling action

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By means of the Fogla Stripping forceps (part number J2892E), hold the peripheral membrane (Fig.6) and start the stripping procedure with a longitudinal movement, making sure not to cause torsions (Fig.7,8). If the adhesion between the membrane and the stroma is too strong, it will cause linear tension on the membrane. It would be advisable at this stage, to perform the stripping in different points, in order to reduce the tension on the tissue, and the endothelial cells mortality.

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


Forceps with tips anatomically designed for the stripping/peeling of the delicate Descemet membrane.

The forceps grasp the membrane at the edge and perform the detachment of the Descemet lenticle by means of a peeling action. At this step the membrane is secured to its support by the suction system.

Cleavage of the membrane is performed with the Fogla hook J2891E

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The membrane must be detached for roughly 90% of the endothelial surface, leaving a small area of peripheral adhesion (Fig.8.1).

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*If the membrane is rolled after the stripping, use a few drops of fluid to spread it out again. Use absorbent pads to remove liquid in excess of the interface between the membrane and the stroma, bringing back the membrane to its original status(fig.9).
Use a dermographic pen to mark the sclera where the membrane is still adherent (fig.10).*

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

Dye the whole tissue using few drops of Trypan Blue for 45 seconds. Remove dye in excess with eye sponges, avoiding any contact with the endothelium. Rinse with PBS drops starting from the adhesion point of the membrane with the stroma. This will avoid the folding of the membrane. Place the tissue in a Petri Dish filled with 2 ml of Sucrose solution 1.8% making sure to firstly lay the side of the sclera marked with the dermographic pen (Fig.11,12). Visually inspect the endothelium.

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

Once the inspection has been carried out, fix a Corneal Claw on the sclera, near the marked area and place the tissue inside the phial as per attached illustrations (Fig.13,14).

This procedure is meant to avoid the lifting of the membrane during the placement inside the phial.

Once placed in the phial, the membrane will position itself in order for the membrane to remain adherent to the stromal layer during transportation (Fig.15).

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Forceps for picking up donor cornea from flask, 1x2 straight, 0.20 mm teeth

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Check the membrane with the slit lamp, if it is detached or partially lifted, repeat the previous procedure. The tissue is now ready for transportation.

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*Procedure edited with the technical and scientific collaboration of the eye bank of Veneto, Director **Diego Ponzin M.D.***

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© Copyright:
E. Janach S.r.l. Via Borgo Vico, 35 - 22100 Como - Italy