

Strumenti chirurgici per oftamologia
e microchirurgia

Ophthalmic Surgical Instruments

ISO 9001 e 13485



E. Janach srl
via Negretti, 3
22100 COMO - Italia
Tel. +39 031 574088
E-mail: janach@janach.it - janach.it

C.F. - P.IVA Reg. Imprese: IT 01557980131
Capitale sociale € 93.600.-int. versati

To whom it may concern

Date: **12/09/2023**

We,

***E. Janach srl
Via Negretti 3
22100, Como, Italy***

declare that our reusable surgical instruments can be sterilized with steam in autoclaves that work with a temperature range of 134°C – 137°C for a minimum of 3 minutes. They can be washed with a thermal disinfection cycle with treatment at least at 90°C for a minimum time of 1 minute.

Marco Nespoli (Legal representative)
E. Janach srl