




The future of vision

# PERFORATED CONFORMER



## Perforated Conformer

|                             |                                                                                                                                                                                                                                                                                                                                                                                                                                        |                             |                            |
|-----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------|----------------------------|
| Brand name                  | AJL PERFORATED CONFORMER                                                                                                                                                                                                                                                                                                                                                                                                               |                             |                            |
| References                  | CF-S-AJL                                                                                                                                                                                                                                                                                                                                                                                                                               | CF-M-AJL                    | CF-L-AJL                   |
| Size                        | Small:<br>22 x 19 x 1.2 mm                                                                                                                                                                                                                                                                                                                                                                                                             | Medium:<br>24 x 21 x 1.2 mm | Large:<br>25 x 22 x 1.2 mm |
| Description and Indications | <p>The AJL PERFORATED CONFORMER is a medical device designed to avoid adherence between the tarsal and bulbar conjunctiva. It is transparent and has two orifices for ventilation at the centre.</p> <p>AJL PERFORATED CONFORMER is used to achieve the separation between tarsal and bulbar conjunctiva. The device prevents the shrinkage of the cavity and maintains the orbital space properly after enucleation/evisceration.</p> |                             |                            |
| Material                    | High resistance acrylic material suitable for ocular use.                                                                                                                                                                                                                                                                                                                                                                              |                             |                            |
| Sterilization method        | Ethylene Oxide (ETO).                                                                                                                                                                                                                                                                                                                                                                                                                  |                             |                            |
| Supply                      | Sterile, single container.                                                                                                                                                                                                                                                                                                                                                                                                             |                             |                            |
| Shelf-life                  | 4.5 years                                                                                                                                                                                                                                                                                                                                                                                                                              |                             |                            |

\* Product with CE marking