


ASPIRA®-aQA

TECHNICAL INFORMATION	
MONOFOCAL MONOBLOC ACRYLIC ASPHERIC 	ASPIRA®-aQA MC 6105 
Type	foldable posterior chamber IOL one-piece with quattro-haptics design
Material	hydrophilic acrylic with UV-absorber; water content 26% at 35°C
Specifics	aspheric anterior surface (aberration free) ⁽²⁾ posterior surface with 360° lens epithelial cell barrier
Diopter Range	
Sphere in 0.5 D steps (optic shape)	10.0 to 30.0 D (biconvex)
Sphere in 1.0 D steps (optic shape)	0.0 to 9.0 D (concave-convex)
Est. A-const ⁽¹⁾	118.7 (optical biometry)

⁽¹⁾ For more details please visit www.humanoptics.com

⁽²⁾ The word „aberration“ as used in this document refers to „spherical aberration“