



PRELOADED SAFELOADER®

INTUITIVE. SIMPLE. FAST.

FOR RELAXED IOL IMPLANTATION

PRECISE AND RELIABLE SUPPLY OF THE IOL IN A UNIQUE AUTOLOADING SYSTEM

- Replaces manual loading of the cartridge by the surgeon or staff
- Intuitive, easy handling
- Reduces the time required for a fast, efficient surgical routine

SEPARATE COMPONENTS FOR MAXIMUM SAFETY

CONTACTLESS ONE-WAY SYSTEM

- Separate, optimal sterilization processes for IOL and injector
- Separate storage of IOL and injector system
- First contact of the IOL with coating additives in the cartridge immediately before implantation

ONE SYSTEM FOR VARIOUS NEEDS

NUMEROUS COMBINATIONS ALLOW INDIVIDUAL FLEXIBILITY RIGHT UP TO THE SURGERY

- Combine the desired premium IOL in the SAFELOADER® with a Mediceal Accuject™ injector of any size¹
- Available for our premium IOLs:



Available lens models	ASPIRA-aA/-aAY	ASPIRA-aXA	DIFF-aA/-aAY	TORICADIFF-aA/-aAY
-----------------------	----------------	------------	--------------	--------------------

Available Accuject™ injector sizes	1.8	2.0	2.2	
------------------------------------	-----	-----	-----	--




INSTANTLY READY TO USE

JUST TRY IT OUT. SIMPLY FINISHED.



Insert injector;
turn to the right



Remove loaded injector;
close cartridge

OUR PREMIUM INTRAOCULAR LENSES

AVAILABLE PRELOADED IN SAFELOADER®



ASPIRA-aXA

**Sub-nano resolution technology**

Premium optic quality

For brilliant, clear, and sharp images

Aspherical optic design

Aberration-free

Suitable for all patients, regardless of corneal aberrations

To preserve depth of field² and enhance contrast sensitivity³**360° LEC barrier**

To prevent PCO

MICS

For minimally invasive implantation

Astigmatism-neutral¹ASPIRA-aA
ASPIRA-aAY **YELLOW**DIFF-aA
DIFF-aAY **YELLOW**TORICADIFF-aA
TORICADIFF-aAY **YELLOW**

HIGH DEFINITION IMPLANTS

PRECISION. OUR PROFESSION.

INNOVATIVE
ENGINEERINGSNR
TECHNOLOGYPREMIUM
QUALITYCUSTOMER
SERVICE¹ Depending on model combination.² Compared to aberration-correcting IOLs.³ Compared to conventional spherical IOLs.