

OASIS

Iris Expander & Collagen Corneal Shield

Expand your view and sustain your access



OASIS IRIS EXPANDER®



A polypropylene ring that gently expands the iris due to insufficiently dilated pupils and sustains visibility throughout the ophthalmic surgery. Each OASIS Iris Expander® is easily inserted and removed. Once removed, the iris returns to its normal shape and function.

- 7mm or 6.25 mm diameter
- Molded polypropylene ring
- Four pockets where Iris is cradled
- Single use, packaged sterile with inserter

Qty	Diameter	Ref#
5 per box	7.0mm	9700
5 per box	6.25mm	9625

SOFT SHIELD® CORNEAL SHIELD



CORNEAL COMFORT & PROTECTION

Provides ocular surface protection following surgery, injection and traumatic or non-traumatic corneal conditions. Available in Quick Shield, 12 hour, 24 hour and 72 hour dissolution times.

- 14mm diameter
- 9mm compound base curve
- Sterile, 10 per box
- Protects the surface of the eye
- Coats the surface of the eye with collagen film as it gradually dissolves
- Certified suitable for use by EDQM
- Hydrates to 70% water content

Description	Ref#
Quick Shield	7002
12 hour Shield	7012
24 hour Shield	7024
72 hour Shield	7072