

Primus-HD®Toric

Primus-HD®Toric Yellow



Preloaded System Hydrophobic Acrylic Toric IOL



Cylinder Power T1 to T6 (1.00 to 4.50 D at IOL Plane)

Toric aspheric 1-piece lens with hydrophobic acrylic material

Glistening free based on in vitro test¹

High stability (median rotation 2°)

100% of the implanted lenses had $\leq 4^\circ$ rotation from end of surgery to 1, 3, 6 months – 1 year postoperatively


1-Step-Preloaded

Model Clear/Yellow	Cylinder Power at IOL Plane	Cylinder Power at Corneal Plane
AT1	1.00 D	0.72 D
AT2	1.50 D	1.07 D
AT3	2.25 D	1.59 D
AT4	3.00 D	2.11 D
AT5	4.00 D	2.80 D
AT6	4.50 D	3.15 D



OPHTHALMO

Pro




Products for ophthalmic surgery

Primus-HD® Toric

Primus-HD® Toric Yellow

Preloaded System / Hydrophobic Acrylic Toric IOL



Technical Data

Model	Primus-HD® Toric Clear/Yellow AT1 – AT6
Optic Design	Posteriorconvex, Toric high order aspheric
Optic & Haptic Materials	Glistening free, hydrophobic Clear with UV-filter, Yellow with blue light filter
Haptic Design	Modified L, 1.5°
Dimension (Optic/OAL)	6.00 mm / 13.00 mm
Power	+5.00 to +36.00 D (in 0.50 D increments)
Cylinder Power	1.00 to 4.50 D (T1 to T6)
Nominal	118.9
Haigis	$a_0 = 1.499$ $a_1 = 0.4$ $a_2 = 0.1$
HofferQ	pACD 5.688
Holladay	sf = 1.949
SRK/T	119.188
SRKII	119.588
Injector	Prosert®, 1-Step-Preload Injector

¹ data on file



Products for ophthalmic surgery

OPHTHALMO Pro GmbH
Im Reihersbruch 1, 66386 Sankt Ingbert
Owner: Thomas Zimmer
Tel.: +49 (68 94) 99 88 770
Fax: +49 (68 94) 99 88 777
Email: office@ophthalmo-pro.de
www.ophthalmo-pro.de/en



www.kestrelphthalmics.co.uk
Email: info@kestrelphthalmics.co.uk
Tel: 01202 696 963