

CROSS LINKING SYSTEM

Available
for daily
hire

Key Specifications

- ∅ Intuitive user friendly interface
- ∅ Light weight, mobile device
- ∅ Suitable for all Cross Linking Protocols
- ∅ 6 – 9 mm beam size
- ∅ UVA intensity adjustment 3-30 mW



The kit is presented
in a secure carry case
for easy everyday
transport

The Intacs XL Cross linking machine is easy to assemble and versatile needing only a small surgical trolley as a base where it can be securely clamped.

CROSS LINKING SINGLE USE PACK^{JCP0007A}

e.janach[®]

The Cross linking single use pack contains all that is needed for a Cross Linking Procedure per patient including:

- Ø Lacrimal cannula 26g
- Ø 5 TNT gauze 10x10cm
- Ø Spatula for epithelium removal
- Ø 2 packs PVA triangular spears
- Ø Vacuum syringe with suction ring, 9mm diameter
- Ø Drape 75x90cm with oval incision film 6x9cm
- Ø Spongystick
- Ø Speculum
- Ø 3ml Syringe
- Ø Bowl 120cc

Also available for purchase separately:

Spatula for epithelial removal (JD2112/10)

Disposable speculum (JD2010/5)

Disposable cxl ring set (JD2112/10)



RIBOFLAVIN SOLUTIONS FOR EPI - OFF & EPI - ON CROSS LINKING PROCEDURES



RIBOFAST Riboflavin

- Ø Dextran Free
- Ø Vitamin E TPGS, corneal penetration enhancer
- Ø EPI- ON, EPI - OFF
- Ø Suitable for all cross linking Protocols



RIBOCROSS te Riboflavin

- Ø Riboflavin solution with 10% Dextran for thinner corneas
- Ø Vitamin E TPGS, corneal penetration enhancer
- Ø EPI- ON, EPI - OFF
- Ø Suitable for all cross linking Protocols

ONLY available for use with Kestrel Ophthalmics UVA Cross Linking Machine

Oasis Collagen Shields available in the following dissolution times:

- Ø Quick Shield (7002)
- Ø 24 hour (7024)
- Ø 12 hour (7012)
- Ø 72 hour (7072)



Contact us

For more information on our products and services please visit our website www.kestrelophthalmics.co.uk or contact us direct at: info@kestrelophthalmics.co.uk

 Kestrel Ophthalmics Ltd

 @KestrelO

Kestrel Ophthalmics Ltd.

Unit 9,
Cabot Business Village,
Holyrood Close, Creekmoor,
Poole, BH17 7BA
01202 696 963